

# Cryogenic Ware

---




**Your Trusted Partner in Life Sciences**



- Abdos cryo vials are manufactured from ultra pure USP class VI certified virgin polypropylene to ensure the highest quality & conforming to US FDA 21 CFR. Free from heavy metals, natural rubber & inhibitory plasticizers.
- Silicon ring meets USP class VI requirement for plastic containers and closures.
- Easy to open & close, the unique grip screw cap provides a secure grip in wet condition.
- Unique combination of plug and flat thread ensuring leak proof and vacuum tight tubes.
- Corrosive and chemical proof white writing area on the tube.
- No leeching with DMSO & Cytotoxin free.
- Suitable for temperatures of -196° to +121° C.
- Specially designed cavity for Cryo Coders for quick & easy visual identification of specimen.
- Star foot design allows for single hand operation when used with Abdos Cryo Box.
- Can withstand pressure upto 95 kpa.
- Sterile tubes has been Validated according to ENISO 11137-2; SAL 10<sup>-6</sup>.

Leeching Performance Report with DMSO




The graph demonstrates no absorbants found to the threshold limits confirming no traceable leeching observed with DMSO.



Code	Product Description	Pack Size	No. Per Case
P60115	ABDOS Cryo Vial External Threaded, 1.0 ml, Star Foot	25 Pcs / Bag	1000
P60116	ABDOS Cryo Vial External Threaded, 1.8 ml, Star Foot	25 Pcs / Bag	1000
P60117	ABDOS Cryo Vial External Threaded, 1.8 ml, Round Bottom	25 Pcs / Bag	1000
P60118	ABDOS Cryo Vial External Threaded, 4.5 ml, Skirted Foot	10 Pcs / Bag	500
P60108	ABDOS Cryo Vial Internal Threaded, 1.0 ml, Star Foot	25 Pcs / Bag	1000
P60109	ABDOS Cryo Vial Internal Threaded, 1.8 ml, Star Foot	25 Pcs / Bag	1000
P60110	ABDOS Cryo Vial Internal Threaded, 1.8 ml, Round Bottom	25 Pcs / Bag	1000
P60112	ABDOS Cryo Vial Internal Threaded, 4.5 ml, Star Foot	10 Pcs / Bag	500
P60113	ABDOS Cryo Vial Internal Threaded, 4.5 ml, Round Bottom	10 Pcs / Bag	500

\* All Abdos Cryo Vials are Gama Radiated Sterile.

\* Abdos range of Cryo Boxes for more details please visit us at [www.abdoslabware.com/cryo-ware.html](http://www.abdoslabware.com/cryo-ware.html)



Cryo Box, PC



Cryo Box System - 100, PC



Freezing Cardboard Cryo Boxes, Cardboard

\* Caution: Do not use Abdos Cryogenic vial in the liquid phase of liquid nitrogen. Improper use may cause liquefied nitrogen to be trapped inside the vial and lead to pressure built up resulting in possible explosion or biohazard release. Can be used in vapour phase of liquid nitrogen.

\* For additional product or technical information, please visit [www.abdoslabware.com](http://www.abdoslabware.com). For distributor information contact [labtech@abdosindia.com](mailto:labtech@abdosindia.com).

## Abdos Labtech Pvt Ltd

720 - 722, DLF Tower A,  
Jasola, New Delhi - 110025, India.

T: +91-11-44800980-986

F: +91-11-44800987

E: [labtech@abdosindia.com](mailto:labtech@abdosindia.com)

### Europe Office

T: +44 - 7958295643

E: [eu.labtech@abdosindia.com](mailto:eu.labtech@abdosindia.com)