

Technical Data Sheet

Doc No–Ab/TDS/UTL-191

Rev No: 02

TECHNICAL DATA SHEET

Abdos Labtech Pvt. Ltd.

Abdos Brand – Skincare Blue Nitrile Gloves – 9.5 inches length

Valid for Catalog No. U20341, U20342, U20343

1. Description/Specification

1.1 Description: Abdos Brand – Skincare Blue Nitrile Gloves – 9.5 inches length

1.1.1 These gloves are ideal product to protect you exposure to bacteria, viruses and chemicals. It also enhances your ability to offer maximum security against cross contamination.

1.1.2 Its unique blue color provides quick, visual differentiation from natural rubber latex and other synthetic gloves.

1.1.3 Ambidextrous, textured fingertips and beaded cuff enhances wet & dry grip.

1.1.4 Custom design enhances comfort & fit.

1.1.5 These gloves are powder free, minimizing the potential for powder related complications.

1.1.6 Conforms to ASTM D6319 and EN 455 Standards.

1.1.7 Skincare Blue Nitrile Gloves Powder free- sample meets the overall migration requirement under Commission Regulation (EU) No 10/2011- Plastic materials and articles shall not transfer their constituents to foodstuff in quantities exceeding 10mg of total constituents released per dm² of food contact surface (mg/dm²) (overall migration limit)

1.2 Design: Hand gloves design.

1.3 Performance Specifications:

| Code | Size | No. Per Case |
|--------|--------|--------------|
| U20341 | Small | 100 |
| U20342 | Medium | 100 |
| U20343 | Large | 100 |



Shiv Ganga Industrial Estate,
Plot no. 1 Lakeshwari, Bhagwanpur,
Roorkee - 247667 Uttrakhand, India
An ISO 9001:2008 Company
A Division of Abdos

Regd. Office :
13A, Government Place (East)
Kolkata – 700069, India
CIN : U85100WB2008PTC123903

P : +91 (0) 1332 235277
E : info@abdoslabware.com
W : www.abdoslabware.com

Technical Data Sheet

1.4 Material: Nitrile gloves

1.5 Color: Blue

1.6 Quality control:


Testing of various characteristics in accordance with the valid Specification

1.7 Other information

1.7.1 These gloves are used for working on lab, machine and any other work.

1.7.2 These gloves are ideal product to protect you any of microbial contamination.

1.7.3 These Flexible gloves are fit in hand & easy to use on working.



Shiv Ganga Industrial Estate,
Plot no. 1 Lakeshwari, Bhagwanpur,
Roorkee - 247667 Uttrakhand, India
An ISO 9001:2008 Company
A Division of Abdos

Regd. Office :
13A, Government Place (East)
Kolkata – 700069, India
CIN : U85100WB2008PTC123903

P : +91 (0) 1332 235277
E : info@abdoslabware.com
W : www.abdoslabware.com