

Technical Data Sheet

Doc No–Ab/TDS/GLP-142

Rev No: 01

TECHNICAL DATA SHEET

Abdos Labtech Pvt. Ltd.

Assorted Workstation Storage Bin, Acrylic

Valid for Catalog No. P70104

1. Description/Specification

1.1 Description: Assorted Workstation Storage Bin, Acrylic.

1.1.1 Assorted Workstation Storage Bin can organize most of their commonly used laboratory disposables in one location.

1.1.2 First section of bin is used to store tubes, vials and pipette tips.

1.1.3 Second section of the bin is used to store and dispense tissue paper, kimwipes and Parafilm M (supplied with super-sharp serrated blade)

1.1.4 Third section of the bins is used to store and dispense gloves.

1.2 Design:

Workstation Storage Bin

1.3 Performance Specifications:

Code	No. of Compartment	L x W x H (mm)	No. Per Case
P70104	5	255 x 135 x 368	1

1.4 Material:

Poly (methyl methacrylate) (PMMA) Acrylic

1.5 Color: Natural

1.6 Quality control:

Testing of various characteristics in accordance with the valid Specification

1.7 Other information:

1.7.1 Strong and lightweight material.



Shiv Ganga Industrial Estate,
Plot no. 1 Lakeshwari, Bhagwanpur,
Roorkee - 247667 Uttarakhand, India
An ISO 9001:2008 Company
A Division of Abdos

Regd. Office :
13A, Government Place (East)
Kolkata – 700069, India
CIN : U85100WB2008PTC123903

P : +91 (0) 1332 235277
E : info@abdoslaware.com
W : www.abdoslaware.com

Technical Data Sheet

1.7.2 Good impact strength, higher than both glass & polystyrene.


1.7.3 Transmits up to 92% of visible light (3 mm thickness)

1.7.4 Better chemical stability.

1.8. Basic Features

1.8.1 Autoclavability: NO, Non autoclavable

1.8.2 Chemical resistance: see website.



Shiv Ganga Industrial Estate,
Plot no. 1 Lakeshwari, Bhagwanpur,
Roorkee - 247667 Uttrakhand, India
An ISO 9001:2008 Company
A Division of Abdos

Regd. Office :
13A, Government Place (East)
Kolkata - 700069, India
CIN : U85100WB2008PTC123903

P : +91 (0) 1332 235277
E : info@abdoslabware.com
W : www.abdoslabware.com