

# Technical Data Sheet

Doc No – Ab/TDS/RB-48

Rev No: 01

## TECHNICAL DATA SHEET

**Abdos Labtech Pvt. Ltd.**

**Weighing Boats, PS**

**Valid for Cat. No. P70501, P70502, P70503, P70504**

### **1. Description/Specification:**

#### **1.1 Description:** Weighing Boats, PS

1.1.1 Weighing Boats are flexible are with smooth surfaces to provide accurate pour-outs with minimal loss.

1.1.2 Biologically inert and resistant to dilute and weak acids, aqueous solution, alcohols & bases.

1.1.3 These boats can withstand temperature between 10°C to 70°C.

1.1.4 Can be used with microwave moisture analyzes.

#### **1.2 Design:**

1.2.1 Weighing Boats are shallow with a wide, flat bottom to resist tipping. They have rounded corners to simplify transfer.

#### **1.3 Performance Specifications:**

##### 1.3.1 Weighing Boats Standard, PS

| Code          | Size   | Colour | L x W x H (mm) | No. per Case |
|---------------|--------|--------|----------------|--------------|
| <b>P70501</b> | Small  | Blue   | 46x46x8        | 250          |
| <b>P70502</b> | Medium | Blue   | 85x85x24       | 250          |

##### 1.3.2 Weighing Boat with Pour Support, PS

| Code          | Size   | Colour | L x W x H (mm) | No. per Case |
|---------------|--------|--------|----------------|--------------|
| <b>P70503</b> | Small  | White  | 43x58x13       | 250          |
| <b>P70504</b> | Medium | White  | 83x132x26      | 250          |



Shiv Ganga Industrial Estate,  
Plot no. 1 Lakeshwari, Bhagwanpur,  
Roorkee - 247667 Uttrakhand, India  
**An ISO 9001:2008 Company**  
A Division of Abdos

**Regd. Office :**  
13A, Government Place (East)  
Kolkata – 700069, India  
**CIN :** U85100WB2008PTC123903

**P :** +91 (0) 1332 235277  
**E :** info@abdoslaware.com  
**W :** www.abdoslaware.com

# Technical Data Sheet

## 1.4 Material:

1.4.2 Complies with requirement of European Pharmacopoeia, 7<sup>th</sup> Edition (2011) & its supplements 7.5 (07/2012) monograph 3.2.2 Plastic container and closures for Pharmaceutical use.

1.4.3 Complies with the requirement of European Pharmacopoeia, 7<sup>th</sup> edition (2011) & its supplements 7.5 (07/2012) monographs 3.2.2 Plastic container and closures for Pharmaceutical use

1.4.4 The product has Passed the USP testing including class VI tests & has been assigned the FDA Drug master file DMF 12123

**1.5 Color:** Blue & White

## 1.6 Quality control:


Testing of various characteristics in accordance with the valid Specification

## 1.7 Basic Features

**1.7.1 Temperature range:** 10<sup>0</sup>C to +70<sup>0</sup>C

**1.7.2 Autoclavable :** No

**1.7.3 Chemical resistance:** See website.



Shiv Ganga Industrial Estate,  
Plot no. 1 Lakeshwari, Bhagwanpur,  
Roorkee - 247667 Uttrakhand, India  
**An ISO 9001:2008 Company**  
A Division of Abdos

**Regd. Office :**  
13A, Government Place (East)  
Kolkata – 700069, India  
**CIN :** U85100WB2008PTC123903

**P :** +91 (0) 1332 235277  
**E :** info@abdoslabware.com  
**W :** www.abdoslabware.com